

# CENTRIFUGE

(digital timer & RPM)

“ISO 9001:2000 CERTIFIED COMPANY”

“CE MARKED PRODUCTS”

Note:- Model no. is connected to the rotor selected. Inter changing of rotor is not permissible as motor are different.

## Model No : CF-611 / CF-622



Capacity	: As per Selected Rotor
Main input	: 220v , 50Hz / 110v , 60Hz
Max. Speed	: As per Selected Rotor
Timer	: Digital Timer (0 to 60 minutes)
Size	: (W) 354 x (D) 454 x (H) 275 mm
Weight	: 14 Kg. Net
Motor	: Universal (AC/DC Carbon Brush Type)
RPM Indicator	: Digital Upto 5 Digits
Speed Regulator	: Step less Electronic

## Model No : CF-611 / RT-03

RT-03



8 x 15ml ANGULAR  
R.P.M : More than 5200  
R.C.F : 3300g  
Material : Aluminum

## Model No : CF-622 / RT-04

RT-04



12 x 15ml ANGULAR  
R.P.M : More than 4900  
R.C.F : 3500g  
Material : Aluminum



LABLINE EQUIPMENTS PVT. LTD.  
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA MACHINERY,  
BARODA 390013, GUJARAT, INDIA.  
PH. NO. +91 265 2632148 / 49  
FAX NO. +91 265 2632151  
E-MAIL : haresh@lablineindia.com

VISIT US AT :- [www.lablineindia.com](http://www.lablineindia.com)