

CENTRIFUGE

Model No	CF - 11	CF - 22
Capacity	4 x 15ml : 4 Nos. Tubes (Swing Out Head)	6 x 15ml : 6 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz 110v,60Hz	: 220v,50Hz / 110v,60Hz
Max. Speed	: 3100 RPM	: 3100 RPM
Max. R.C.F.	: 1750 x g	: 1750 x g
W x D x H	: 300 x 300 x 245mm	: 300 x 300 x 245mm
Weight	: 3.1 Kg. Net	: 3.2 Kg. Net
Motor	: Universal (AC/DC Carbon Brush type)	: Universal (AC/DC Carbon Brush type)



Without TIMER



Auto Stop Model

Model No	CF - 33	CF - 44
Capacity	4 x 15ml : 4 Nos. Tubes (Swing Out Head)	6 x 15ml : 6 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz / 110v,60Hz	: 220v,50Hz / 110v,60Hz
Max. Speed	: 3100 RPM	: 3100 RPM
Max. R.C.F.	: 1750 x g	: 1750 x g
Timer	: 0 to 30 min.	: 0 to 30 min.
W x D x H	: 300 x 300 x 245mm	: 300 x 300 x 245mm
Weight	: 3.4 Kg. Net	: 3.5 Kg. Net
Motor	: Universal (AC/DC Carbon Brush type)	: Universal (AC/DC Carbon Brush type)

CENTRIFUGE (ABS BODY With Locking LID)

Model No	CF - 501	CF - 502
Capacity	4 x 15ml : 4 Nos. Tubes (Swing Out Head)	6 x 15ml : 6 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz 110v,60Hz	: 220v,50Hz 110v,60Hz
Max. Speed	: 3200 RPM	: 3200 RPM
Max. R.C.F.	: 1600 x g	: 1700 x g
W x D x H	: 315 x 315 x 268mm	: 315 x 315 x 268mm
Weight	: 3.9 Kg. Net	: 4.0 Kg. Net
Motor	: Universal (AC/DC Carbon Brush type)	: Universal (AC/DC Carbon Brush type)
Body	: ABS Body (Light Wt.) With Metal Shielding	: ABS Body (Light Wt.) With Metal Shielding



Without TIMER