

CENTRIFUGE (ABS BODY With Locking LID)

Most reliable unit for routine centrifuging tests of medical practitioners, pathological laboratories and hospitals provided with 5 (Five) step speed.

Model No : CF - 503

CF - 504

	8 x 15ml	12 x 5ml
Capacity	: 8 Nos. Tube (Swing Out Head)	: 12 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz 110v,60Hz	: 220v,50Hz 110v,60Hz
Max. Speed	: 3200 RPM	: 3200 RPM
Max. R.C.F.	: 1800 x g	: 1400 x g
W x D x H	: 315 x 315 x 268mm	: 315 x 315 x 268mm
Weight	: 4.1 Kg. Net	: 4.1 Kg. Net
Motor	: Universal (AC/DC Carbon Brush type)	: Universal (AC/DC Carbon Brush type)
Body	: ABS Body (Light Wt.) With Metal Shielding	: ABS Body (Light Wt.) With Metal Shielding



Without TIMER



Auto Stop Model

Model No : CF - 506

CF - 507

	4 x 15ml	6 x 15ml
Capacity	: 4 Nos. Tubes (Swing Out Head)	: 6 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz 110v,60Hz	: 220v,50Hz 110v,60Hz
Max. Speed	: 3200 RPM	: 3200 RPM
Max. R.C.F.	: 1600 x g	: 1700 x g
W x D x H	: 315 x 315 x 268mm	: 315 x 315 x 268mm
Weight	: 4.0 Kg. Net	: 4.1 Kg. Net
Motor	: Universal (AC/DC Carbon Brush type)	: Universal (AC/DC Carbon Brush type)
Timer	: 0 to 30 Minutes	: 0 to 30 Minutes
Body	: ABS Body (Light Wt.) With Metal Shielding	: ABS Body (Light Wt.) With Metal Shielding

Quality is not an accident, it is an intelligent efforts