

CENTRIFUGE (ABS BODY With Locking LID)

Model No : CF - 508

CF - 509

	8 x 15ml	12 x 5ml
Capacity	: 8 Nos. Tube (Swing Out Head)	: 12 Nos. Tubes (Swing Out Head)
Mains Input	: 220v,50Hz 110v,60Hz	: 220v,50Hz 110v,60Hz
Max. Speed	: 3200 RPM	: 3200 RPM
Max. R.C.F.	: 1800 x g	: 1400 x g
W x D x H	: 315 x 315 x 268mm	: 315 x 315 x 268mm
Weight	: 4.2 Kg. Net	: 4.2 Kg. Net
Motor	: Universal (AC/DC Carbon Brush type)	: Universal (AC/DC Carbon Brush type)
Timer	: 0 to 30 Minutes	: 0 to 30 Minutes
Body	: ABS Body (Light Wt.) With Metal Shielding	: ABS Body (Light Wt.) With Metal Shielding



Auto Stop Model

HEMATOCRIT CENTRIFUGE (With Digital Timer & Digital RPM)

NEW



Model No : HCF-12 & HCF-24

Capacity	: 12 / 24 Nos. Capillary Tube
Mains Input	: 220v,50Hz / 110v,60Hz
Max. Speed	: 12000 RPM (Digital)
Timer	: 0 to 60 Min. (Digital)
W x D x H	: 354 x 454 x 230mm
Weight	: 12.0 Kg. Net
Motor	: Universal (AC/DC Carbon Brush type)
Lock	: Electrical Lock
RPM Indicator	: Digital Upto 5 Digits
Speed Regulator	: Step Less

The Customer is not a hurdle in business, He is purpose of Business