




"ISO 9001:2008 & EN ISO 13485:2012"  
CERTIFIED COMPANY  
"CE MARKED PRODUCTS"

## DRY BATH INCUBATOR

[ FIX TEMPERATURE TYPE ]

### Model No : DBI - 22

|                  |                                     |
|------------------|-------------------------------------|
| Mains Input      | : 230V/50Hz   110V/60Hz A.C. Supply |
| W x D x H        | : 195 (W) x 86 (D) x 80 (H)         |
| Body             | : ABS & Metal Body                  |
| Weight           | : 1.5 Kg. Net                       |
| Temp. Controller | : Digital Type                      |
| Temp. Fix        | : 37°C                              |
| Accuracy         | : $\pm 0.5^{\circ}\text{C}$         |
| Heating Block    | : 16 Holes of 12.7mm Diameter       |




## DRY BATH INCUBATOR

### Model No : DBI - 11

|                  |                                                           |
|------------------|-----------------------------------------------------------|
| Mains Input      | : 230V/50Hz   110V/60Hz A.C. Supply                       |
| W x D x H        | : 168 (W) x 232 (D) x 93 (H)                              |
| Body             | : Metal Body                                              |
| Weight           | : 2.3 Kg. Net                                             |
| Temp. Controller | : Digital Type                                            |
| Temp. Range      | : 37°C To 70°C                                            |
| Accuracy         | : $\pm 0.5^{\circ}\text{C}$ @ 37°C, $\pm 3\%$ Except 37°C |
| Heating Block    | : 24 Holes of 12.7mm Diameter<br>or As Selected - any one |
| Timer            | : 1 To 999 ( Digital Type )                               |




Block - 1




12.7 MM 24 HOLE

Block - 2




10.5 MM 28 HOLE

Block - 3



18.5 MM - 4 HOLE  
12.7 MM - 6 HOLE  
15.5 MM - 5 HOLE  
10.5 MM - 7 HOLE

Block - 4



12.7 X 12.7 MM  
20 SQUAR HOLE