

NEEDLE DESTROYER (ND-66)

Model No : ND - 66



Mains Input	: 230v, 50Hz 110v, 60Hz
Body	: Metal Body
Gauge Of Needles	: Can Burn 30 to 20 Gauge
Destruction Time	: 1.5 to 3 Seconds
Burning Temp.	: 1400° C to 1800° C
Weight	: 1.9 Kg. Net
Overall Size	: 155 x 128 x 80 mm

Finest quality with incredible price

- ☞ Premium quality transformer with CRNGO lamination and copper winding [same as ND 22].
- ☞ S.S. Hard and grind syringe cutter [same as ND 22].
- ☞ CRCA steel body with finest quality powder coating.
- ☞ Most popular brand name of “ LABLINE ”.

All this for long lasting safety

You must have services support facility for minor problems although Technical Guidance, Drawing, Diagram, Spares and further Technical Assistance will be provided by us.



LABLINE EQUIPMENTS PVT. LTD.
614 / 1 GIDC, MAKARPURA, B/H. HIMALAYA M ACHINERY,
BARODA 390010, GUJARAT, INDIA.
PH. NO. +91 265 3927448 / 49 / 50
FAX NO. +91 265 2632151
E-MAIL : haresh@lablineindia.com

VISIT US AT :- www.lablineindia.com