


MOTOR BRUSH LESS WITH AC DRIVE



Model :- CF-888



Model :- PRP - 01

- ➔ ROTOR :- 8 TUBE SWING OUT
- ➔ POSITION :- 90 DEGREE HORIZONTAL TYPE
- ➔ TUBE SIZE :- 1 - DIA 13 x 75 MM LONG - 5ML / 2 - DIA 13 x 100 MM LONG - 7ML
 3 - DIA 16 x 100 MM LONG - 15 ML / 4 - DIA 16 x 120 MM LONG - 15ML
- ➔ R.C.F :- 2000 g
- ➔ RPM :- 3500 max



Model :- SC-101



Model :- PRP - 02


- ➔ ROTOR :- 4 TUBE SWING OUT
- ➔ POSITION :- 90 DEGREE HORIZONTAL TYPE
- ➔ TUBE SIZE :- 1 - DIA 13 x 75 MM LONG - 5ML / 2 - DIA 13 x 100 MM LONG - 7ML
 3 - DIA 16 x 100 MM LONG - 15 ML / 4 - DIA 16 x 120 MM LONG - 15ML
- ➔ R.C.F :- 2350 g
- ➔ RPM :- 4000 max

MOTOR UNIVERSAL (AC/DC CARBON BRUSH TYPE)



Model :- CF-600

- ➔ ROTOR :- 8 TUBE SWING OUT
- ➔ POSITION :- 90 DEGREE HORIZONTAL TYPE
- ➔ TUBE SIZE :- 1 - DIA 13 x 75 MM LONG - 5ML / 2 - DIA 13 x 100 MM LONG - 7ML
3 - DIA 16 x 100 MM LONG - 15 ML / 4 - DIA 16 x 120 MM LONG - 15ML



Model :- PRP - 01

- ➔ R.C.F :- 2000 g
- ➔ RPM :- 3500 max



Model :- CF-513

- ➔ ROTOR :- 4 TUBE SWING OUT
- ➔ POSITION :- 90 DEGREE HORIZONTAL TYPE
- ➔ TUBE SIZE :- 1 - DIA 13 x 75 MM LONG - 5ML / 2 - DIA 13 x 100 MM LONG - 7ML
3 - DIA 16 x 100 MM LONG - 15 ML / 4 - DIA 16 x 120 MM LONG - 15ML



Model :- PRP - 02

- ➔ R.C.F :- 2350 g
- ➔ RPM :- 4000 max